

CWB-001

Column Flushing System




CWB-001

Column Flushing System

● Product Overview

- » Cleaning and maintenance of the column after use, no longer need to occupy the liquid chromatography instrument alone to flush the column, and improve the efficiency of liquid chromatography.
- » Column Pre-balancing: pre-balancing of the column is carried out before the experiment to save balancing time during the experiment.
- » Equipped with a 5-channel online automatic switching valve, it can automatically realize the flushing or pre-balancing of 5 columns.
- » Optional 4-channel cleaning solvents, freely switch solvents, and easily change gradients to meet the needs of flushing or pre-balancing of different columns.
- » Full-color intelligent touch screen control, or controlled by PC, powerful software support, can realize automatic switching of different flushing methods (switching flow rate, ratio, solvent, etc.), automatic shutdown, etc.

● Featured Product Recommendations:



ELSD



HPLC



上海通微分析技术有限公司
Unimicro (Shanghai) Technologies Co., Ltd.

上海通微分析技术有限公司

电话: 021-38953588 38953390 38953570

邮箱: info@unimicrotech.com.cn

地址: 上海市浦东新区张江高科技园区松涛路489号
C01座

网址: www.unimicrotech.com.cn



Unimicro Technologies

美国通微技术股份有限公司

Unimicro Technologies, Inc.

Tel.: 925-846-8638 Fax: 925-401-9548

Address: 23046 Avenida de la Carlota, Suite 600,
Laguna Hills CA 92653

Website: www.unimicrotech.com